

CLICK CALL CORD



Although call cords may be the simplest type of nurse call device, they may be too small for patients with limited hand dexterity; and if patients aren't sure they pressed the call button, they won't know if they called the nurse.

Curbell Medical's Click Call Cord addresses these issues. It's inexpensive, UL 1069-listed, and designed to improve the patient experience.

The call button makes an audible and tactile click, so patients can feel assured that they called the nurse. It fits naturally in the hand, with a raised lip to help prevent the device from slipping out of a patient's hand – making it ideal for the elderly and those with weak hand strength.

- ▶ Audible and tactile call button
- ▶ Raised lip to prevent device from slipping out of a patient's hand
- ▶ UL 1069 listed



Plug Type	Nurse Call System	Length	Part Number	Part Number
1/4"	Universal	6 inches	CCE-0000-000	29395
1/4"	Universal	8 feet	CCE-0000-800	29190
1/4"	Universal	10 feet	CCE-0000-A00	29191
1/4"	Universal	12 feet	CCE-0000-C00	29192
1/4"	Universal	20 feet	CCE-0000-H00	29396

©2026 Curbell Medical Products, Inc. MAP1834C