

# ASSISTIVE CALL DEVICES

Our Assistive Call Devices offer ease-of-use and comfort to those with special needs. They feature ultra-sensitive touch switches, making them ideal for patients with limited manual dexterity, and are compatible with most nurse call systems.



## ZEPHYR™

- ▶ Breath-activated device for patients with compromised manual dexterity
- ▶ Single (Sip or Puff), and Dual (Sip-and-Puff) models available
- ▶ The reusable base (gooseneck, clamp, and cables) is easily cleaned between patients
- ▶ The disposable switch assembly is changed between patients to help with infection control
- ▶ Arm can be adjusted to almost any position (adjustable clamp fits side rail or headboard)



## BREATHCALL®

- ▶ Breath-activated device for patients with compromised manual dexterity
- ▶ Single and dual models available
- ▶ Adjustable clamp fits side rail or headboard
- ▶ Arm can be adjusted to almost any position
- ▶ Additional straws and filters available
- ▶ Single patient use



## E-Z CALL™

- ▶ Can be placed anywhere for activation including the cheek, foot, hand, etc.
- ▶ Large tented design makes it easy for patient to see
- ▶ Linen clip securely attaches to bedding
- ▶ Velcro strip included for side rail mounting



## MECHANICAL PAD

- ▶ Perfect for patients who may have difficulty using standard call cords or pillow speakers
- ▶ Large sensitive surface area
- ▶ Activates from any point on the pad's surface
- ▶ Linen clip securely attaches to bedding



## PADCALL®

- ▶ Large call surface (3" diameter)
- ▶ Sensitive, air-activated switch
- ▶ Gripper feet hold pendant in place
- ▶ Linen clip securely attaches to bedding
- ▶ Custom lengths available



## PRESSCALL®

- ▶ Air activated
- ▶ Fits comfortably in patient's palm
- ▶ Available in dual and single styles
- ▶ Linen clip securely attaches to bedding
- ▶ Custom lengths available