




MAGNETIC BED BREAKAWAY CABLE

In busy healthcare environments, beds may be moved before staff remember to unplug them, potentially damaging the bed connection, cable, wall receptacle, or patient station. Curbell's Magnetic Breakaway Bed Cable maintains a secure, reliable connection, but separates safely to help minimize this damage.



The unique design is quick to connect and easy to align, with hidden magnets that ensure a proper connection while reducing the risk of attracting objects that could cause a short circuit. The cable's recessed, pinless connector has no moving parts, unlike other magnetic breakaway cables with spring-loaded “pogo” pins which frequently break or jam. Its blue and orange housing makes it easy to identify among the many cables often found behind the bed.



BED END	LENGTH	WALL END	LENGTH	ANGLE	MODEL	PART
 37-Pin female, Hillrom*	6 ft.	 37-Pin male	1 ft.	straight	BCC-3802-A01-A06-001	30780
				90°	BCC-3803-A01-A06-001	30784
	1 ft.	 37-Pin male	6 ft.	straight	BCC-3802-A06-A01-001	30782
				90°	BCC-3803-A06-A01-001	30786
 37-Pin male, Stryker**	6 ft.	 37-Pin male	1 ft.	straight	BCC-3802-A01-E06-001	30781
				90°	BCC-3803-A01-E06-001	30785
	1 ft.	 37-Pin male	6 ft.	straight	BCC-3802-A06-E01-001	30783
				90°	BCC-3803-A06-E01-001	30787

* Also for Linet®, Joerns®, and Novum® beds

** Also for Amico®, Arjo®, Carroll® (GHC®), Kreg Medical, Next Health®, OnCare®, Rotec, Sizewise®, Stryker® London, Ufiano Medical, and VitalGo beds

©2026 Curbell Medical Products, Inc. MAP1962A

BREAKAWAY CABLE LOCKING SOLUTIONS


For a cable to separate at the breakaway point, it needs to be secured at both the bed end and the wall end.



WPA-0719
Inline Cable Wall Plate/Bracket
Allows you to secure inline breakaway cables.



WPA-0713
90° Cable Wall Plate/Bracket
Allows you to secure 90° breakaway cables.



WPA-0608
Bracketless Wall Plate
Allows you to secure the breakaway cable using the bed cable mounting screws.



BRK-1007-001
Dual 90° Locking Bracket
Allows you to secure two 90° breakaway cables to an existing dual 37-pin wall plate.